

KEY FEATURES	440	330	220	110	BENEFIT
Performance	xxxxxx	xxxxx	xxxx	xx	
Processing channels	15	12	10	6	For accurate processing and a more precise sound experience and speech intelligibility
Fine-tuning channels	15	12	10	6	More options for personalised sound to cater for user preferences and hearing needs
WidexLink wireless connectivity	●	●	●	●	Fast, echo-free transmission of sound between hearing aids and DEX devices
SoundSense Learn	●	●	●		Machine learning algorithm to effortlessly guide the user to a better, more personalised hearing experience
SoundSense Adapt	●	●	●		Adaptive learning to remember the user's preferred settings across multiple parameters for each environment, resulting in better sound
Fluid Sound Controller	●	●	●	●	Adjustment of more parameters in every sound class and in real time for more natural sound
Fluid Sound Analyzer (sound classes)	11 (IE)	7 (IE)	4	3	Intuitive analysis of sound environment in order to prioritise and select the most important sound sources
Programs	5	4	3	3	Wide range of listening programs for specific and challenging situations
Smartwind Manager	●				Speech and other sounds audible in wind
High-frequency boost	●				Clear and crisp sound in high frequencies (in music, for example)
Speech Enhancer RT	RT/IE	IE			Reduces noise and enhances speech in real time in response to individual sound situations and hearing loss
Digital Pinna	●	●			Helps locate and focus naturally on the voice in front or behind by mimicking the normal pinna effect
HD Locator	●	●	●		For better speech intelligibility in noise with channel-specific directionality for speech emphasis
TruSound Softener	●	●	●		Preserves comfort for sudden sounds at any input level by taking the edge off impact/impulse sounds
Soft-level noise reduction	●	●	●	●	Unwanted soft background noises are reduced for greater listening comfort
Noise Reduction	●	●	●	●	Listening comfort in noise and more speech
ZEN IE	●	●	●	●	A tinnitus management tool for reducing the effects of tinnitus using soothing tones and/or noise; additional programs with Zen+
Audibility Extender	●	●	●	●	Improves perception and production of important high-frequency speech sounds
Variable Speed Compression	●	●	●	●	Automatic gearing for smooth sound and stable audibility
EVOKE App	●	●	●	●	Easy connection and control from any smartphone (F2 model)
TONELINK App	●	●	●	●	Easy connection and control from any smartphone (all hearing aid models except CIC-M)